

SOME THINGS ARE A
TOUCH BETTER.
SEMPERCARE PREMIUM.



sempercare[®]  **premium**

EXAMINATION GLOVES BY SEMPERMED

SEMPERCARE PREMIUM

EXTRA COMFORT MADE BY SEMPERMED.

For your professional top performance the best working material is just about good enough. Of course this credo is also true for our examination gloves. With the Sempercure Premium we have developed a product which meets the highest expectations and sets new standards.

- For the Sempercure Premium glove a new extra wide porcelain dipping mould with longer fingers was developed. The optimized length of the fingers significantly improves the **wearing comfort** of our powder-free natural latex examination glove.
- Our powder-free natural latex examination glove excels by its extraordinary, **well-gliding** surface. The textured fingers provide a good grip which is indispensable for **safe instrument handling**.
- The **synthetic inner lining** achieved by the special Sempermed technology **makes it easy and fast to don**.
- Like all products made by Sempermed, the Sempercure Premium convinces with its particular **skin friendliness**.

TECHNICAL DATA



sempercure®  premium

TYPE	unsterile examination glove for single use, powder-free
MATERIAL	natural latex (NR)
COLOUR	natural white
GLOVE SHAPE	flat with rolled rim
OVERALL LENGTH/SIZE as per EN 455-2	XS, S, M, L, XL 240 mm
WALL THICKNESS measured in double layer in the palm area	min. 0.16 mm
IMPERMEABILITY/FORCE AT BREAK as per EN 455-1/as per EN 455-2	AQL 1.5 / $\geq 9N$
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years
PACKAGING	unsterile, in dispenser carton: XS, S, M, L 100 pcs./XL 90 pcs. in transport carton: XS, S, M, L 1.000 pcs./XL 900 pcs./PU
MARKING	EN 1041, EN 980 and ISO 15223
PURPOSE as set out in MDD 93/42 EEC	Medical device class I

CONTACT US!

Semperit Technische Produkte GmbH · Division Sempermed
Austria · 1031 Vienna · Modecenterstraße 22 · Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
e-mail: sempermed@semperit.at · www.sempermed.com