

**POWDER-
FREE**

FOR SENSITIVE HANDS

The Sempermed Derma PF offers good reasons to switch to powder-free gloves.

HIGHLY SKIN-FRIENDLY

This powder-free glove provides improved safety for staff, patients and visitors as there are no powder-related dermatological, allergic and toxic health risks.

EASY TO DON

With its synthetic inner coating the glove is easy to put on due to its patented net structure.

SURE GRIP

The Sempermed Derma PF has more grip and thus combines safety, tactical properties and instrument guidance.

GOOD VALUE FOR MONEY

The long-term cost reduction is significant, if you take the very high follow-up costs of powder damage into consideration. The cost of purchasing Sempermed Derma PF are minimal by comparison.



sempermed®  **DERMA PF**

DESIGNED AROUND YOUR HANDS.

TECHNICAL DATA



sempermed® DERMA PF

TYPE	Sterile surgical glove, powder-free with synthetic lining		
MATERIAL	Natural rubber latex (NR)		
COLOUR	White		
GLOVE SHAPE	Fully anatomical with rolled rim		
OVERALL LENGTH/SIZE as per EN 455-2	5 ½, 6 and 6 ½	265 mm	
	7, 7 ½ and 8	275 mm	
	8 ½ and 9	280 mm	
WALL THICKNESS measured in single layer in the palm area	0.16 - 0.24 mm		
IMPERMEABILITY as per EN 455-1	AQL 1.5		
FORCE AT BREAK as per EN 455-2	≥ 12 N		
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	min. 3 years		
STERILISATION	Gamma radiation with at least 2.5 Mrad (25 kGy)		
PACKAGING	left and right glove with turned up cuff in fibre-free internal pouch, ozone-tight, sealed in medical peel pack.		
	in dispenser carton with sterilisation indicator: 50 pairs		
	in transport carton with sterilisation indicator: 300 pairs		
ARTICLE NUMBERS	Size 5.5	822751522	Size 7.5 822751722
	Size 6.0	822751602	Size 8.0 822751802
	Size 6.5	822751622	Size 8.5 822751822
	Size 7.0	822751702	Size 9.0 822751902
MARKING	EN 1041, EN 980, ISO 15223		

CONTACT US!

Semperit Technische Produkte GmbH · Division Sempermed

1031 Vienna · Modocenterstraße 22 · Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630

e-mail: sempermed@semperit.at · www.sempermed.com