

# FOR PARTICULARLY SENSITIVE AREAS

## Light Absorbing Colour

The unique brown colour of this powder-free Latex surgical glove reduces light reflection during applications in the operating theatre. The lower reflection levels help reduce the strain on eyes during sensitive procedures.

## Ideal Tactile Properties

The extra-thin surgical glove Sempermed Senso allows for precise fine motor activities as well as ideal tactile sense. In addition to the ideal tactile sense Sempermed Senso has an even wall thickness, which feels like an almost intangible second skin.



## A Sure Grip

A perfect surgical glove must never be slippery. At the same time, instruments must not "stick to the fingers". With its micro-rough surface texture the Sempermed Senso guarantees a good grip – even in tricky situations.

## Maximum Comfort

Developed specifically for the requirements of microsurgery, its fully anatomical form with slightly curved fingers corresponds with the natural hand position. This allows for fatigueless work and a calm, relaxed handling of instruments – even during long operations.

**sempermed**<sup>®</sup> 

**senso**

DESIGNED AROUND YOUR HANDS.

# TECHNICAL DATA



**sempermed®**  **senso**

<b>TYPE</b>	Sterile surgical glove, powder-free with synthetic lining		
<b>MATERIAL</b>	Natural latex (NR)		
<b>COLOUR</b>	Brown		
<b>GLOVE SHAPE</b>	Fully anatomical with rolled rim		
<b>OVERALL LENGTH/SIZE</b> as per EN 455-2	5 ½	250 mm	
	6, 6 ½, 7, 7 ½		
	8, 8 ½ and 9	280 mm	
<b>WALL THICKNESS</b> in palm area	0.14 - 0.18 mm		
<b>IMPERMEABILITY</b> as per EN 455-1	AQL 1.0		
<b>FORCE AT BREAK</b> as per EN 455-2	≥ 12 N		
<b>DURABILITY</b> in original packaging if stored as per DIN 7716, ISO 2230	min. 3 years		
<b>STERILISATION</b>	Gamma radiation with at least 2.5 Mrad (25 kGy)		
<b>PACKAGING</b>	Left and right hand with turned up cuff in fibre-free internal pouch, ozone-tight, sealed in medical peel pack. In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs		
<b>ARTICLE NUMBERS</b>	Size 5.5	821751521	Size 7.5 821751721
	Size 6.0	821751601	Size 8.0 821751801
	Size 6.5	821751621	Size 8.5 821751821
	Size 7.0	821751701	Size 9.0 821751901
<b>MARKING</b>	EN 1041, EN 980, ISO 15223		

## CONTACT US!

Semperit Industrial Products Ltd.  
25 Cottesbrooke Park · Heartlands  
Daventry · Northants · NN11 8YL  
Tel.: +44 1327 313140 · Fax: +44 1327 313149  
sales@semperit.co.uk · www.sempermed.com

Semperit Technische Produkte Gesellschaft m.b.H.  
Division Sempermed  
1031 Wien · Modecenterstraße 22  
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630  
sempermed@semperit.at · www.sempermed.com