

NITRILE SKIN²

DISCOVER YOUR 2ND SKIN

PROTECTION FROM INFECTIONS – virus-proof for users and patients
EXCELLENT SENSE OF TOUCH – through optimised wall thickness
HIGH WEARING COMFORT – through the innovative, soft material

Protection against infections
Excellent sense of touch

High wearing comfort
Free from natural latex

sempercare[®]  **nitrile skin²**

EXAMINATION GLOVES BY SEMPERMED

YOUR NEEDS UNDER CONTROL

In addition to protection, wearing comfort and skin-friendliness, a good sense of touch is the most important criterion when choosing the right examination glove. The new Sempercure Nitrile skin² combines these properties in a unique, exceptionally light glove.

The soft material developed by our European Research Centre allows for maximum sense of touch and comfort due to optimised wall thickness, and at the same time maximum protection from infections.

Discover your 2nd skin!

- **Protects** users and patients from the transmission of infections. In addition to the impermeability requirements of the European standard, the gloves also fulfil the virus impermeability requirements of ASTM F 1671.
- With the **newly developed material**, the Sempercure Nitrile Skin² feels extremely **soft and light**.
- The elastic material offers highest wearing comfort and a maximum **sense of touch**.
- The structured fingertips guarantee a **sure grip** at all times, even in a moist environment.
- The glove is powder-free and does **not** contain any **natural latex proteins**. This makes it ideal for people with type I allergy who do not want to miss the comfort of latex.
- The tried and tested **Sempercure washing process** during production of the glove guarantees maximum skin-friendliness of the Sempercure Nitrile Skin², both for the patients and for the users.
- The **violet colour** makes it easy to distinguish from natural latex gloves and allows for easy recognition.

TECHNICAL DATA



sempercure®  nitrile skin²

TYPE	Unsterile, powder-free examination glove for single use
MATERIAL	Nitrile butadiene rubber (NBR)
COLOUR	Violet
GLOVE SHAPE	Flat with rolled rim
OVERALL LENGTH/SIZE as per EN 455-2	XS, S, M, L, XL 240 mm
WALL THICKNESS measured in double layer in the palm area	min. 0.10 mm
IMPERMEABILITY / FORCE AT BREAK as per EN 455-1 / as per EN 455-2	AQL 1.5 / 6 N
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years
PACKAGING	In dispenser carton XS, S, M, L 200 gloves loose, XL 180 gloves loose 10 dispenser cartons per shipping carton
MARKING	EN 1041, EN 980, EN 420, ISO 15223
PURPOSE as set out in MDD 93/42 EEC as set out in PPE 89/686 EEC	Medical device class I Single-use protective glove for minimal risks

CONTACT US!

Semperit Technische Produkte GmbH · Division Sempermed
Austria · 1031 Vienna · Modecenterstraße 22 · Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
e-mail: sempermed@semperit.at · www.sempermed.com