

**MORE
GRIP**

A NEW GRIP TO YOUR WORLD

The new Sempermed Supreme+ provides you with more grip in every situation while combining maximum safety, ideal tactical properties and perfect instrument guidance.

ABSOLUTELY SKIN-FRIENDLY

The specially developed, automatic finishing process contains intensive and multiple washing steps. Skin-irritating substances resulting from the production process and the natural latex protein concentration are thus reduced to the technically lowest level.

THE LINING SYSTEM

The special net-like lining structure optimises the contact surface between the film and the skin, thus considerably reducing the friction resistance when putting on the glove. The result: The new Sempermed Supreme+ facilitates easy donning of the gloves.

DOUBLE DONNING

For your safety and easy donning use Sempermed Supreme as under glove and Sempermed Supreme+ as outer glove.

Recommendation: If you wear different sizes, put on the larger pair first. Or use gloves of the same size.



sempermed®  **supreme+**

DESIGNED AROUND YOUR HANDS.

TECHNICAL DATA



sempermed® supreme+

TYPE	Sterile surgical glove, powder-free, with synthetic lining			
MATERIAL	Natural latex (NR)			
COLOUR	White			
GLOVE SHAPE	Fully anatomical with rolled rim			
OVERALL LENGTH/SIZE as per EN 455-2	5 ½, 6 and 6 ½	270 mm	7, 7 ½ and 8	280 mm
	8 ½ and 9	285 mm		
WALL THICKNESS measured in single layer in the palm area	0.19 - 0.24 mm			
IMPERMEABILITY as per EN 455-1	AQL 1.0			
FORCE AT BREAK as per EN 455-2	≥ 12 N			
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	min. 3 years			
STERILISATION	Gamma radiation with at least 2.5 Mrad (25 kGy)			
PACKAGING	Left and right hand with turned up cuff in fibre-free internal pouch, ozone-tight, sealed in medical peel pack. In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs			
ARTICLE NUMBER	Size 5.5	822851521	Size 7.5	822851721
	Size 6.0	822851601	Size 8.0	822851801
	Size 6.5	822851621	Size 8.5	822851821
	Size 7.0	822851701	Size 9.0	822851901
MARKING	EN 1041, EN 980, ISO 15223			

CONTACT US!

Semperit Industrial Products Ltd.
25 Cottesbrooke Park, Heartlands · Daventry Northants NN11 8YL
Tel.: +44 1327 313140 · Telefax: +44 1327 313148
sales@semperit.co.uk · www.sempermed.com

Semperit Technische Produkte GmbH · Division Sempermed · 1031 Wien · Modecenterstraße 22 · www.sempermed.com