

sempermed® 
classic

CLASSIC

THE ESTABLISHED CLASSIC

IDEAL TACTILE PROPERTIES – highly developed dipping process
HIGH SAFETY – ideal compounding
EXCELLENT GRIP – high comfort



THE ESTABLISHED CLASSIC

■ IDEAL TACTILE PROPERTIES

The highly developed dipping process guarantees constant wall thickness from fingers to cuff. This results in ideal tactile sensitivity and less hand fatigue during even longer operations.

■ SKIN CARE

The glove is lightly powdered with pure cornstarch which makes it easy to don. The multiple washing process (leaching) makes the natural latex glove skin friendly and very comfortable to wear.

■ HIGHER SAFETY

Sempermed's manufacturing process complies with all international quality standards and guarantees the required Acceptable Quality Level of 1.5. The compounding of the Natural Rubber Latex and production technology result in a higher force at break and therefore in a higher safety.

■ EXCELLENT GRIP AND COMFORT

Sensitivity, good control and a safe grip are essential for successful surgery. This glove has a micro-rough non-slip surface to enhance its ability to grip.



TECHNICAL DATA

sempermed® classic

TYPE	Sterile surgical glove for single use, lightly powdered			
MATERIAL	Natural rubber latex (NR)			
COLOUR	White			
GLOVE SHAPE	Fully anatomical with rolled rim			
OVERALL LENGTH/SIZE as per EN 455-2	6 and 6 ½		min. 260 mm	
	7, 7 ½ and 8		min. 270 mm	
	8 ½ and 9		min. 280 mm	
WALL THICKNESS measured in single layer in the palm area	min. 0.16 mm			
IMPERMEABILITY / FORCE AT BREAK as per EN 455-1 / as per EN 455-2	AQL 1.5 / ≥12 N			
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	5 years			
STERILISATION	Gamma radiation with at least 2.5 Mrad (25 kGy)			
PACKAGING	<p>1. Sterile: Left and right glove with turned up cuff in fibre-free inner pouch, sealed in medical peel pack. In dispenser carton with sterilisation indicator: 70 pairs In transport carton with sterilisation indicator: 420 pairs</p> <p>2. Bulk non sterile: 500 left hand and 500 right hand gloves loose. Each 500 packed in a PE-bag in one transport carton.</p>			
ARTICLE NUMBERS	Size 6.0	826054600	Size 8.0	826054800
	Size 6.5	826054620	Size 8.5	826054820
	Size 7.0	826054700	Size 9.0	826054900
	Size 7.5	826054720		
MARKING	EN 1041, EN 980, ISO 15223			

CONTACT US!

Semperit Technische Produkte Gesellschaft m. b. H. Division Sempermed
Austria · 1030 Vienna · Modecenterstraße 22 · Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
e-mail: sempermed@semperit.at · www.sempermed.com