

# SENSO FOR PARTICULARLY SENSITIVE AREAS

IDEAL TACTILE PROPERTIES – extra-thin wall thickness  
BROWN COLOUR – reduced strain on eyes  
A SURE GRIP – optimal instrument guidance



# FOR PARTICULARLY SENSITIVE AREAS

## ■ LIGHT ABSORBING COLOUR

The unique brown colour of this powder-free latex surgical glove reduces light reflection during applications in the operating theatre. The lower reflection levels help reduce the strain on eyes during sensitive procedures.

## ■ IDEAL TACTILE PROPERTIES

The extra-thin surgical glove Sempermed Senso allows for precise fine motor activities as well as ideal tactile sense. In addition to the ideal tactile sense Sempermed Senso has an even wall thickness, which feels like an almost intangible second skin.

## ■ A SURE GRIP

A perfect surgical glove must never be slippery. At the same time, instruments must not "stick to the fingers". With its micro-rough surface texture the Sempermed Senso guarantees a good grip – even in tricky situations.

## ■ MAXIMUM COMFORT

Developed specifically for the requirements of microsurgery, its fully anatomical form with slightly curved fingers corresponds with the natural hand position. This allows for fatigueless work and a calm, relaxed handling of instruments – even during long operations.



## TECHNICAL DATA

### sempermed® senso

<b>TYPE</b>	Sterile surgical glove for single use, powder-free			
<b>MATERIAL</b>	Natural rubber latex (NR)			
<b>COLOUR</b>	Brown			
<b>GLOVE SHAPE</b>	Fully anatomical with rolled rim			
<b>OVERALL LENGTH/SIZE</b> per EN 455-2	5 ½		min. 250 mm	
	6, 6 ½, 7, 7 ½, 8, 8 ½ and 9		min. 280 mm	
<b>WALL THICKNESS</b> in palm area	0.14 - 0.18 mm			
<b>IMPERMEABILITY / FORCE AT BREAK</b> as per EN 455-1 / as per EN 455-2	AQL 1.0 / ≥12 N			
<b>DURABILITY</b> in original packaging if stored as per DIN 7716, ISO 2230	3 years			
<b>STERILISATION</b>	Gamma radiation with at least 2.5 Mrad (25 kGy)			
<b>PACKAGING</b>	Left and right hand with turned up cuff in fibre-free inner pouch, ozone-tight, sealed in medical peel pack.			
	In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs			
<b>ARTICLE NUMBERS</b>	Size 5.5	821751521	Size 7.5	821751721
	Size 6.0	821751601	Size 8.0	821751801
	Size 6.5	821751621	Size 8.5	821751821
	Size 7.0	821751701	Size 9.0	821751901
<b>MARKING</b>	EN 1041, EN 980, ISO 15223			

### CONTACT US!

Semperit Technische Produkte Gesellschaft m. b. H. Division Sempermed  
Austria · 1030 Vienna · Modecenterstraße 22 · Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630  
e-mail: [sempermed@semperit.at](mailto:sempermed@semperit.at) · [www.sempermed.com](http://www.sempermed.com)