



A NEW GENERATION OF SEMPERMED SURGICAL GLOVES

SEMPERMED GREEN M

with its unique colour reduces the reflections of light during applications in the operating theatre. This new generation of gloves thus is perfectly designed for sensitive procedures which need minimisation of eye-strain as well as a good contrast to organs.

A sure grip

The ideal combination of strength-to-thickness ratio allows precise fine motor activities as well as ideal tactile sense. A micro-roughened surface ensures extra sensitivity in the fingertips for optimum safety in perfect instrument handling, even in a moist environment. This makes the glove an excellent choice for use in ophthalmic, cardiac, micro-vascular and many other delicate operations.

Double Donning

The synthetic lining of the Sempermed GREEN M gloves makes them easy to don even with moist hands. In particularly difficult situations or with high risk patients, two gloves worn on top of each other provide maximum safety. Recommendation: If you wear different sizes, put on the larger pair first or use gloves of the same size.

sempermed®  **GREEN M**

DESIGNED AROUND YOUR HANDS.

TECHNICAL DATA



sempermed® GREEN M

TYPE	Sterile surgical glove, powder-free with synthetic lining	
MATERIAL	Natural latex (NR)	
COLOUR	Green	
GLOVE SHAPE	Fully anatomical with rolled rim	
OVERALL LENGTH/SIZE as per EN 455-2	5 ½ 6, 6 ½, and 7 7 ½, 8, 8 ½ and 9	275 mm 280 mm 290 mm
WALL THICKNESS in palm area	min. 0.2 mm	
IMPERMEABILITY as per EN 455-1	AQL 1.0	
FORCE AT BREAK as per EN 455-2	≥ 12 N	
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	min. 3 years	
STERILISATION	Gamma radiation with at least 2.5 Mrad (25 kGy)	
PACKAGING	Left and right hand with turned up cuff in fibre-free internal pouch, ozone-tight, sealed in medical peel pack. In dispenser carton with sterilisation indicator: 40 pairs In transport carton with sterilisation indicator: 240 pairs	
ARTICLE NUMBERS	Size 5.5 823756521 Size 6.0 823756601 Size 6.5 823756621 Size 7.0 823756701	Size 7.5 823756721 Size 8.0 823756801 Size 8.5 823756821 Size 9.0 823756901
MARKING	EN 1041, EN 980, ISO 15223	

CONTACT US!

Semperit Technische Produkte Gesellschaft m.b.H. · Division Sempermed
1031 Wien · Modecenterstraße 22 · Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
e-mail: sempermed@semperit.at · www.sempermed.com