



Higher puncture resistance and better chemical resistance

For users with type I allergies

More comfort

The blue colour makes it easier to distinguish from latex gloves.

Perfect fit and sensitivity

**sempercare<sup>®</sup>**  **nitrile**

DESIGNED AROUND YOUR HANDS.

## SEMPERCARE NITRILE

# THE SYNTHETIC LATEX ALTERNATIVE

The tried and tested safety of a nitrile examination glove – combined with wearing properties similar to latex – make the Sempercare Nitrile the new generation for an optimum glove solution without natural latex.

- The new Sempercare Nitrile is even softer and more elastic, thus offering a **perfect fit** and **sensitivity**.
- The soft material makes it more **comfortable** to wear and safer to work with, even for longer periods.
- The glove is powder-free and does not contain any natural latex proteins. This makes it perfect for users with **type I allergies** who do not want to miss the greater comfort of elastic natural latex.
- The multiple washing process during production makes the glove particularly **skin friendly**.
- The structured fingertips guarantee a **good grip**, even in moist conditions
- The higher **puncture resistance** and improved **chemical resistance** of nitrile – with the same thickness as natural latex – offer the user protection while maintaining comfort.
- The blue colour makes it **easier to distinguish from latex gloves**.



## TECHNICAL DATA

### sempercare® nitrile

<b>TYPE</b>	Powder-free examination glove for single use
<b>MATERIAL</b>	Nitrile butadiene rubber (NBR)
<b>COLOUR</b>	blue
<b>GLOVE SHAPE</b>	flat with rolled rim
<b>SIZE / OVERALL LENGTH</b> as per EN 455-2	S, M, L, XL 240 mm
<b>WALL THICKNESS</b> measured in double layer in the palm area	min. 0.16 mm
<b>IMPERMEABILITY / FORCE AT BREAK</b> as per EN 455-1 / as per EN 455-2	AQL 1.5 / $\geq 6$ N
<b>DURABILITY</b> in original packaging if stored as per DIN 7716, ISO 2230	3 years
<b>PACKAGING</b>	1. unsterile: in dispenser carton S, M, L 100 gloves loose, XL 90 gloves S, M, L 200 gloves loose, XL 180 gloves in transport carton: S, M, L 1.000 gloves loose, XL 900 gloves S, M, L 2.000 gloves loose, XL 1800 gloves 2. sterile: S, M, L 50 pairs in dispenser carton, gloves with turned up cuff in fibre-free internal pouch, sealed in medical peel pack in transport carton: S, M, L 300 pairs/PU
<b>MARKING</b>	EN 1041, EN 980, ISO 15223, EN 420 (unsterile)
<b>PURPOSE</b> as set out in MDD 93/42 EEC and PPE 89/686 EEC	Medical device class I (unsterile + sterile) Single-use protective glove for minimal risks (unsterile)

#### FOR SALES SUPPORT AND INFORMATION PLEASE CONTACT:

##### Semperit Industrial Products Ltd.

25 Cottesbrooke Park · Heartlands · Daventry · Northants · NN11 8YL

Tel: +44 1327 313140 · Fax: +44 1327 313149 · sales@semperit.co.uk · www.sempermed.com

Semperit Technische Produkte GmbH · Division Sempermed · A-1031 Vienna · Modecenterstraße 22 · www.sempermed.com