

S semper
med

velvet XP



**Fits like a
second skin**



Hospital



Dental



Primary
care



Long-term
care



Veterinary

CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
9 Straits View, #08-10A Marina One West Tower
Singapore 018937, Singapore
Tel.: +65 6274 4861 · Fax: +65 6274 6977

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348

S semper
med

velvet XP



**Fits like a
second skin**



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
9 Straits View, #08-10A Marina One West Tower
Singapore 018937, Singapore
Tel.: +65 6274 4861 · Fax: +65 6274 6977

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348



velvet XP



**Fits like a
second skin**



Hospital



Dental



**Primary
care**



**Long-term
care**



Veterinary

CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
9 Straits View, #08-10A Marina One West Tower
Singapore 018937, Singapore
Tel.: +65 6274 4861 · Fax: +65 6274 6977

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348



In addition to protection, wearing comfort, and skin-friendliness, a good sense of touch is another important criteria when choosing the right examination glove. The sempermed® velvet XP combines these factors and is a trusted solution for a wide range of applications.

- **Protects:** users and patients from the transmission of infections according to EN 455-1. In addition to the impermeability requirements of the European standard, the gloves also fulfil the virus impermeability requirements of ASTM F 1671.
- The elastic material provides **high wearing comfort** and tactile sensitivity.
- The structured fingertips guarantee a **secure grip**, even in a moist environment.
- The glove is powder-free and does **not** contain any **natural latex proteins**.
- The sempermed® velvet XP is both labelled as medical device (class 1) and single-use protective glove (PPE category III). Therefore, the glove is suitable for use in areas such as laboratory/ pharma, food/catering, professional cleaning/hygiene.
- The **violet color** makes it easy to distinguish from natural latex gloves and allows for easy recognition.

Technical data

Type	non-sterile, powder-free examination glove for single use
Color	violet blue
Size / Overall length as per EN 455-2	XS, S, M, L, XL / median 240 mm
Wall thickness measured in single layer in the palm area	min. 0.06 mm typical 0.07 mm
Barrier performance as per EN 455-1	AQL 1.0
Material	Nitrile Butadiene Rubber (NBR)
Glove surface	finger textured
Glove shape	flat with rolled rim
Force at break as per EN 455-2	median ≥ 6 N typical 7 N
Shelf life in original packaging if stored according to product specification	5 years
Labeling & conformity to standards	EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 2859-1, ASTM D6319, ASTM F1671/F1671M

Purpose

Single-use medical device class I as set out in MD Regulation (EU) 2017/745

Single-use protective glove – PPE category III (protection against chemical substances for limited time)* as set out in PPE Regulation (EU) 2016/425

Suitable for food contact as set out in Regulation (EC) 1935/2004 on Food Contact Materials



Packaging

Dispenser box XS, S, M, L 100 / XL 90 pieces
Dispenser box XS, S, M, L 200 / XL 180 pieces
10 dispenser boxes per shipping carton

Article numbers

XS	100pc	3000015971	200pc	3000015991
S	100pc	3000015972	200pc	3000015992
M	100pc	3000015973	200pc	3000015993
L	100pc	3000015974	200pc	3000015994
XL	90pc	3000015975	180pc	3000015995

DM1/12 Version 1 10/2024 EN