

New Look & Specifications

Upgraded Sterile Examination Glove Packaging

We are proud to present the refreshed packaging for our sterile examination gloves. The medical examination gloves previously known as *sempercure*, are now under the *sempermed* brand, ensuring a consistent look and feel across our entire glove family while maintaining the same trusted quality.

		<p>Medical examination gloves now under the sempermed umbrella brand, replacing sempercure brand</p> <p>NEW</p>
		<p>Packaging design to incorporate our new branding</p> <p>NEW</p> <p>Our new progressive logo inspired by the infinity symbol, reflects our long-term focus on quality, reliability, & customer service. Semper [sem.per] = latin; always</p>
<p>Semperit Investments Asia Ptd. Ltd.</p>	<p>HARPS Investment Asia Pte. Ltd.</p>	<p>Manufacturer name and address</p> <p>NEW</p>
<p>Semperit Technische Produkte GmbH</p>	<p>EU REP according to MDR (EU) 2017/745 HARPS Europe GmbH</p>	<p>EU Representative, also serving as the Regulatory Importer for EU markets</p> <p>NEW</p>
<p>LOT + YYYY-MM XXXXXX YYYY-MM</p>	<p>UDI (xx) xxxxxxxxxxxxxx (xx) xxxxxx (xx) xxxxxx</p>	<p>Implementation of data matrix printing according to MDR UDI on all packaging levels</p> <p>SAME</p>
	<p>For detailed product-specific information, please refer to the following page(s).</p>	<p>Article numbers and barcodes</p> <p>NEW</p>
<p>AQL ≤ 1.5</p>	<p>AQL AQL ≤ 1.0</p>	<p>Change of AQL for pinholes in the sterile examination glove portfolio</p> <p>UPDATE</p> <p>(please refer to individual product specifications for more information)</p>
		<p>UDI now visible on two sides for improved logistics</p> <p>UPDATE</p>
<p>The product contains dithiocarbamates – not to be used in the case of hypersensitivity to these substances.</p>	<p>Allergen IV</p>	<p>Now featuring the new EN 455-3:2023 symbol (replacing the previous warning text for easier recognition and consistent understanding)</p> <p>UPDATE</p>
		<p>Same conformity to regulations and product classification / same Basic UDI-DIs (please refer to individual product specifications and Declarations of Conformity)</p>

sempermed® safe sterile

Great things now come in sempermed packaging!



Our strong sterile nitrile examination gloves are now branded as **sempermed® safe sterile** and will be exclusively available with an AQL 1.0 for pinholes.

Please refer to the following reference numbers, GTIN codes:

Reference number	Article Description	GTIN code Dispenser	GTIN code Carton	Packing	Dispenser Dimensions	Carton Dimensions	Load ability
sempermed® safe sterile <i>previously sempercare nitrile sterile</i>							
3000016225	sempermed® safe sterile S	9001570551818	9001570551801	50 pairs/ dispenser (S, M, L, XL), 300 pairs per carton	L 155 mm x W 135 mm x H 250 mm	L 480 mm x W 285 mm x H 261 mm	20ft container: 840 cartons
3000016226	sempermed® safe sterile M	9001570551658	9001570551641				40ft container: 1690 cartons
3000016227	sempermed® safe sterile L	9001570551672	9001570551665				40ft HC container: 1920 cartons
3000016228	sempermed® safe sterile XL	9001570551696	9001570551689				



Hospital



Dental



Primary
care



Long-term
care



Veterinary



Please contact your Area Sales Manager in case of any questions!

DM3/9_19 Version 1 10/2025 EN

find more
information online



CONTACT US!
sempermed@harpsglobal.com
www.sempermed.com

HARPS Investment Asia Pte. Ltd.
9 Straits View, #08-10A Marina One West Tower
Singapore 018937, Singapore
Tel.: +65 6274 4861 · Fax: +65 6274 6977

HARPS Europe GmbH
Wiedner Guertel 9-13,
1100 Vienna, Austria
Tel.: +43 1 9346348