

EDITION IC THE ALL-ROUNDER

POLYMER COATED – for easy donning

TEXTURED – for strong grip

HIGH WEARING COMFORT – due to natural latex material

Powder-free latex examination glove

Polymer coated – for easy donning

Textured – for good handling of instruments

excellent elasticity

THE ALL-ROUNDER

The Sempermed® Edition IC combines a micro rough surface on the outside with a polymer on the inside of the glove. With its excellent grip properties of the micro rough surface, the Sempermed® Edition IC is ideal for handling instruments for instance in the lab as well as for normal daily nursing tasks. The special polymer inner coating makes the latex glove very comfortable to wear and easy to don.

- With special synthetic inner coating
- The surface ensures good gripping and secure instrument handling
- The high elasticity of natural latex makes the glove very comfortable to wear

TECHNICAL DATA

TYPE	Powder free examination and disposable protective glove, non sterile	
MATERIAL	Natural latex (NR)	
COLOUR	Natural white	
GLOVE SHAPE	Flat with rolled rim, ambidextrous	
SIZE/OVERALL LENGTH as per EN 455-2	XS, S, M, L, XL / median 240 mm	
WALL THICKNESS measured in double layer in the palm area	Min. 0.20 mm Typical 0.26 mm	
BARRIER PERFORMANCE as per EN 455-1	AQL 1.5	
FORCE AT BREAK as per EN 455-2	Median \geq 6 Newton Typical production value \geq 7 Newton	
DURABILITY in original packaging if stored as per ISO 2230	3 years	
PACKAGING	In dispenser carton:	XS,S,M,L 100 pcs. / XL 90 pcs.
	In transport carton:	XS,S,M,L 1.000 pcs. / XL 900 pcs.
LABELLING	EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5, EN 455-1, EN 455-2, EN 455-3, EN 455-4, ISO 15223-1, ISO 2859-1, ASTM D3578 (except stress at 500% elongation), ASTM F1671	
PURPOSE as set out in MD Regulation (EU) 2017/745 as set out in PPE Regulation (EU) 2016/425 as per Reg. (EC) 1935/2004 on Food Contact Materials	Medical device class I Single-use protective glove – PPE category III (protection against chemical substances for limited time)* Suitable for food contact	



CONTACT US!

Semperit Investments Asia Pte Ltd
8 Jurong Town Hall Road,
#29-03 to 06 JTC Summit, Singapore 609434
✉ sempermed@semperitgroup.com
🌐 www.sempermed.com

Semperit Technische Produkte Gesellschaft m.b.H.
Division Sempermed
Austria · 1030 Vienna · Modecenterstrasse 22
☎ +43 1 79 777-0 📠 +43 1 79 777-630
✉ sempermed@semperitgroup.com 🌐 www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2019 Sempermed