

# SAFE+ EXTRA LONG. EXTRA SAFE.

ADDED SPLASH PROTECTION – due to long cuff  
BARRIER PROTECTION\* – good Chemical Resistance  
CONFIDENT GRIP – thanks to Finger Texturing

30 cm length provides added splash protection  
Finger texturing offers a confident grip  
Protects against certain chemicals\*  
Nitrile material is suitable for sensitive skin  
Striking blue colour

## EXTRA LONG. EXTRA SAFE.

The Sempercure® Safe+ effectively safeguards your hands when added barrier as well as splash protection are required. Thanks to its 30cm-long cuff, wrists and forearms are better protected, which makes it a suitable choice for glove users who work with risk patients or handle substances that can be harmful to the skin.

- **ENHANCED SAFETY:**  
additional splash protection thanks to the extra-long cuff
- **CHEMICAL RESISTANCE:**  
barrier protection against certain liquids and chemicals
- **WORK WITH CONFIDENCE:**  
finger texturing and good tactile sensitivity
- **SUITABLE FOR SENSITIVE SKIN:**  
nitrile material helps to prevent type-1 allergies

## TECHNICAL DATA

<b>TYPE</b>	powder-free disposable examination glove	
<b>MATERIAL</b>	nitrile-Butadien-Rubber (NBR)	
<b>COLOUR</b>	blue	
<b>GLOVE SHAPE</b>	flat with rolled rim	
<b>SIZE/OVERALL LENGTH</b> as per EN 455-2	XS, S, M, L, XL / median 300 mm	
<b>WALL THICKNESS</b> measured in double layer in the palm area	0.20 mm	
<b>BARRIER PERFORMANCE</b> Tightness as per EN 455-1	AQL 1.5	
<b>FORCE AT BREAK</b> as per EN 455-2	≥ 8 N	
<b>DURABILITY</b> in original packaging if stored as per DIN 7716, ISO 2230	3 years	
<b>PACKAGING</b>	in dispenser box: XS, S, M, L 100 gloves loose, XL 90 pieces in transport box: XS, S, M, L 1000 gloves loose, XL 900 pieces	
<b>LABELLING</b>	EN 455-1/-2/-3/-4, EN 420, EN 421, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5, ISO 15223-1	
<b>PURPOSE</b> as set out in MD Regulation (EU) 2017/745 as set out in PPE Regulation (EU) 2016/425	Medical device class I Single-use protective glove – PPE category III (protection against chemical substances for limited time)*	
<b>ARTICLE NUMBERS</b>	Size XS    8267 58031	Size S    8267 58033
	Size M    8267 58035	Size L    8267 58037
	Size XL    8267 58039	



## CONTACT US!

Semperit Investments Asia Pte. Ltd,  
8 Jurong Town Hall Road, #29-03 to 06 The JTC Summit,  
Singapore 609434, Singapore  
Tel.: +65 6274 4861 · Fax.: +65 6274 6977  
sempermed@semperitgroup.com · www.sempermed.com

Semperit Technische Produkte Gesellschaft m.b.H.  
Triester Bundesstraße 26, 2632 Wimpassing, Austria  
Tel.: +43-2630-310-0  
sempermed@semperitgroup.com · www.sempermed.com

\*Please note the Sempermed Chemical Resistance list is available on [www.sempermed.com/sempermed\\_informs/chemical\\_resistance](http://www.sempermed.com/sempermed_informs/chemical_resistance)

**IMPORTANT NOTE:** The latest product information is available at [www.sempermed.com](http://www.sempermed.com) under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2022 Sempermed