

VELVET

SOFT AS VELVET AND LIGHT AS A FEATHER

YOUR DAILY CHOICE – a thin & light glove with a force at break of $\geq 6\text{ N}^*$
COMFORTABLE FIT – thanks to the soft glove design
ENHANCED SENSE OF TOUCH – good finger dexterity

High wearing comfort
Soft glove material

Enhanced sense
of touch



SEMPERCARE® VELVET

SOFT AS VELVET AND LIGHT AS A FEATHER

The comfortable and light Sempercure® Velvet is a soft latex-free option for professionals looking for high wearing comfort and an enhanced sense of touch combined with effective protection against infections.

- Feather light and velvety soft, the Sempercure® Velvet offers good wearing comfort and tactile sensitivity.
- The classification as medical device (class I) as set out in as set out in MD Regulation (EU) 2017/745 and as protective glove (category III) as set out in PPE Regulation (EU) 2016/425 makes the Sempercure® Velvet suitable for a variety of applications.
- Sempercure® Velvet is powder free (as per EN 455-3) and does not contain any natural rubber protein. As such, this glove is perfectly suitable for glove users with Type-I allergy.
- Individual size color coding on dispenser boxes allows for the easy identification of the correct glove size.
- Packed in 100 or 200 pieces per dispenser box.



TECHNICAL DATA

TYPE Unsterile, powder-free examination glove for single use

MATERIAL Nitrile Butadiene Rubber (NBR)

COLOUR Lavender blue

GLOVE SHAPE Flat with rolled rim

SIZE/OVERALL LENGTH XS, S, M, L, XL / median 240 mm
as per EN 455-2

WALL THICKNESS min. 0.10 mm
measured in double layer in the palm area typical 0.11-0.12 mm

BARRIER PERFORMANCE AQL 1.5
as per EN 455-1

FORCE AT BREAK ≥ 6 Newton
as per EN 455-2

DURABILITY 3 years
in original packaging if stored
as per DIN 7716, ISO 2230

PACKAGING Available in 100 and 200 pieces per dispenser box
In dispenser box XS, S, M, L 100 / XL 90 pieces
In dispenser box XS, S, M, L 200 / XL 180 pieces
10 dispenser boxes per shipping carton

LABELLING EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4,
EN ISO 374-5, EN 455 1-4, EN ISO 15223-1

PURPOSE Medical device class IIa
Single-use protective glove – PPE category III
(protection against chemical substances for limited time)*
as set out in MD Regulation (EU) 2017/745
as set out in PPE Regulation (EU) 2016/425



ARTIKELNUMMERN		100 pcs.		Size M		Size XL	
		Size XS	300 000 6795	300 000 6797	300 000 6799		
	200 pcs.	Size S	300 000 6796	Size L	300 000 6798		
		Size XS	300 000 6800	Size M	300 000 6802	Size XL	300 000 6804
		Size S	300 000 6801	Size L	300 000 6803		

CONTACT US!

Semperit Investments Asia Pte Ltd
8 Jurong Town Hall Road,
#29-03 to 06 JTC Summit, Singapore 609434
✉ sempermed@semperitgroup.com
🌐 www.sempermed.com

Semperit Technische Produkte Gesellschaft m.b.H.
Division Sempermed
Austria · 1030 Vienna · Modecenterstrasse 22
☎ +43 1 79 777-0 📠 +43 1 79 777-630
✉ sempermed@semperitgroup.com 🌐 www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2020 Sempermed

