

SEMPERCLEAN MC QUALITY YOU CAN FEEL

EASY TO DON – as a result of the synthetic inner coating

SECURE GRIP – thanks to micro-rough surface and outer surface treatment

HIGH WEARING COMFORT – due to anatomical shape

QUALITY YOU CAN FEEL

The Semperclean® MC is a specialized protective glove made of natural latex and is suitable to the needs of medical-technical areas, such as laboratories, the pharmaceutical industry and R&D. The synthetic lining that makes the Semperclean® MC so easy to don is the result of many years of development and the web-like structure of the glove surface reduces friction resistance, improving the wearing comfort considerably. The fully anatomical shape with curved fingers helps to prevent hand fatigue and facilitates optimized instrument handling.

KEY FEATURES:

- Synthetic inner coating for easy donning
- A secure grip due to special outer surface treatment
- Anatomical shape reduces hand fatigue
- Available in 6 sizes, including intermediate sizes for optimal fit
- Recommended usage: laboratory, pharmaceutical industry, research & development



TECHNICAL DATA

TYPE	disposable protective and examination glove
MATERIAL	natural rubber latex (NR)
COLOUR	natural white
GLOVE SHAPE	fully anatomical with rolled rim
SIZE/OVERALL LENGTH as per EN 455-2	5½, 6 and 6½ / median 270 mm 7, 7½ and 8 / median 280 mm 8½ and 9 / median 285 mm
WALL THICKNESS measured in double layer in the palm area	0,42 mm ± 0,05 mm
BARRIER PERFORMANCE as per EN 455-1, EN 374-2	AQL 1.5
FORCE AT BREAK as per EN 455-2 typical production values	median 6 N median 12 N
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years
PACKAGING	unsterile in dispenser carton: 50 pairs in transport carton: 500 pairs
LABELLING AND CONFORMITY WITH STANDARDS	EN 420 / EN ISO 21420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5, EN 455 1-4, EN ISO 15223-1
PURPOSE as set out in MD Regulation (EU) 2017/745 as set out in PPE Regulation (EU) 2016/425	medical device CE class I single-use protective glove – PPE CE category III (protection against chemical substances for limited time)*



CONTACT US!

Semperit Investments Asia Pte Ltd
8 Jurong Town Hall Road,
#29-03 to 06 JTC Summit, Singapore 609434
✉ sempermed@semperitgroup.com
🌐 www.sempermed.com

Semperit Technische Produkte Gesellschaft m.b.H.
Division Sempermed
Austria · 1030 Vienna · Modecenterstrasse 22
☎ +43 1 79 777-0 📠 +43 1 79 777-630
✉ sempermed@semperitgroup.com 🌐 www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2021 Sempermed