

sempermed® 
derma plus

DERMA PLUS PRODUCED SAFETY

HIGH TEAR RESISTANCE – through double-immersion production process

PERFECT GRIP – through micro-rough surface

ANATOMICALLY PERFECT SHAPE – with broader back of the hand
and naturally curved finger



SEMPERMED® DERMA PLUS

PRODUCED SAFETY

Manufactured while using the specially developed double-dipping process the Sempermed® Derma plus excels with good donning properties and an AQL of 0.65.

PRODUCT HIGHLIGHTS

- High tear resistance and reliable impermeability through double-immersion production process.
- Good grip through micro-rough surface texture, even in a moist environment.
- Anatomically perfect shape with broader and more comfortable back of the hand and naturally curved fingers.



TECHNICAL DATA

sempermed® derma plus

TYPE	Sterile surgical glove for single use, lightly powdered		
MATERIAL	Natural rubber latex (NR)		
COLOUR	Natural white		
GLOVE SHAPE	Fully anatomical with rolled rim		
SIZE/OVERALL LENGTH as per EN 455-2	5 ½, 6 and 6 ½ 7, 7 ½ and 8 8 ½ and 9	median 270 mm median 280 mm median 285 mm	
WALL THICKNESS measured in single layer in the palm area	0.19 – 0.24 mm		
BARRIER PERFORMANCE Tightness exceeding EN 455-1 (AQL 0.65 instead of AQL 1.5)	AQL 0.65		
FORCE AT BREAK exceeding EN 455-2 (≥ 12 N instead of 9 N)	≥ 12 N		
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	5 years		
STERILISATION	Radiation STERILE R		
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in medical peel pack. In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs		
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223		
PURPOSE according to MP Regulation 93/42/EWG according to PPE Regulation (EU) 2016/425	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time)*		
ARTICLE NUMBERS	Size 5 ½ 823152529 Size 6 823152609 Size 6 ½ 823152629 Size 7 823152709	Size 7 ½ 823152729 Size 8 823152809 Size 8 ½ 823152829 Size 9 823152909	



CONTACT US!

Semperit Investments Asia Pte Ltd
8 Jurong Town Hall Road,
#29-03 to 06 JTC Summit, Singapore 609434
E-Mail: sempermed@semperitgroup.com
www.sempermed.com

Semperit Technische Produkte Gesellschaft m. b. H.
Division Sempermed
Austria · 1030 Vienna · Modecenterstrasse 22
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
E-Mail: sempermed@semperitgroup.com · www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/sempermed_informs/chemical_resistance

IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the 'Products' tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2019 Sempermed